

Cleaning Procedure for CarryOn™



To clean the CarryOn™ use only a damp cloth with a little mild detergent and dry thoroughly afterwards with a clean, soft cloth/towel.

DO NOT USE:

- Bleach
- Solvents (acetone, white spirit, mineral turpentine, ethanol (methylated spirits), Xylene etc.)
- Abrasive cleaners
- Scouring pads, steel wool or other abrasive material
- Dripping-wet cloths – to protect from electrical hazard.

Use of these substances/methods will void the warranty.